

Anti-Bax Rabbit antibody



Catalog Number: 167534

Product name

Anti-Bax Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bax Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Bax

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of Bax in C6,3T3,Hela cell lysates using Bax Rabbit pAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.