

Anti-FUBP1/FBP Rabbit antibody



Catalog Number: 167537

Product name

Anti-FUBP1/FBP Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

FUBP1/FBP Rabbit polyclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human FUBP1

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal



Western blot detection of FUBP1 in Jurkat cell lysates using FUBP1 Rabbit pAb(1:1000 diluted). Predicted band size:68KDa. Observed band size:74KDa.