

Anti-Phospho-Rb (Ser780) Rabbit antibody



Catalog Number: 167543

Product name

Anti-Phospho-Rb (Ser780) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, IP

Antibody description

Phospho-Rb (Ser780) Rabbit polyclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser780 of human Rb

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Validations



Western blot detection of Phospho-Rb (Ser780) in K562,C6 cell lysates using Phospho-Rb (Ser780) Rabbit pAb(1:1000 diluted).Predicted band size:106kDa.Observed band size:110kDa.