

# Anti-KAT2A/GCN5 Rabbit antibody



Catalog Number: 167547

## Product name

Anti-KAT2A/GCN5 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

KAT2A/GCN5 Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

IP: 1/20


## Preparation

Antigen: Recombinant protein of human GCN5

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Western blot detection of KAT2A / GCN5 in K562 cell lysates using KAT2A / GCN5 Rabbit pAb(1:1000 diluted).Predicted band size:94kDa.Observed band size:94kDa.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal