

# Anti-TSG101 Rabbit antibody



Catalog Number: 167554

## Product name

Anti-TSG101 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TSG101 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human TSG101

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/2000

IP: 1/20

## Validations



Western blot detection of TSG101 in K562 cell lysates using TSG101 Rabbit pAb(1:1000 diluted).Predicted band size:44KDa.Observed band size:44KDa.