

# Anti-USP10 Rabbit antibody

Catalog Number: 167559

## Product name

Anti-USP10 Rabbit antibody

## Specificity

Human

## Antibody description

USP10 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human USP10

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of USP10 in HeLa cell lysates using USP10 Rabbit pAb(1:1000 diluted). Predicted band size:87kDa. Observed band size:110kDa.



Immunocytochemistry of USP10(green) in HeLa cells using USP10 Rabbit pAb at dilution 1/50, and DAPI(blue)