

Anti-UQCRB Rabbit antibody



Catalog Number: 167560

Product name

Anti-UQCRB Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

UQCRB Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human UQCRB

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/50

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of UQCRB in Rat Brain lysates using UQCRB Rabbit pAb(1:1000 diluted). Predicted band size:14kDa. Observed band size:14kDa.