

# Anti-Ataxin 3 Rabbit antibody



Catalog Number: 167568

## Product name

Anti-Ataxin 3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Ataxin 3 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Ataxin 3

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Ataxin 3 in HeLa, A549, HL-60, U251, U87-MG cell lysates using Ataxin 3 Rabbit mAb(1:1000 diluted). Predicted band size: 42kDa. Observed band size: 42kDa.