

Anti-ATF3 Rabbit antibody



Catalog Number: 167574

Product name

Anti-ATF3 Rabbit antibody

Specificity

Human, Mouse

Antibody description

ATF3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human ATF3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of ATF3 in MCF-7,U251,HL-60 cell lysates using ATF3 Rabbit mAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.