

Anti-Tissue Factor Rabbit antibody



Catalog Number: 167576

Product name

Anti-Tissue Factor Rabbit antibody

Specificity

Human

Antibody description

Tissue Factor Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Tissue Factor

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of Tissue Factor in A549,HL-60,U2OS cell lysates using Tissue Factor Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:45kDa.