

Anti-TFEB Rabbit antibody



Catalog Number: 167580

Product name

Anti-TFEB Rabbit antibody

Specificity

Human

Antibody description

TFEB Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human TFEB

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of TFEB in A549,HL-60 cell lysates using TFEB Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:65-70kDa.