

Anti-MELK Rabbit antibody



Catalog Number: 167589

Product name

Anti-MELK Rabbit antibody

Specificity

Human

Antibody description

MELK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MELK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of MELK in U2OS,C6 cell lysates using MELK Rabbit mAb(1:1000 diluted).Predicted band size:75kDa.Observed band size:75kDa.