

# Anti-Ly-6A/E Rabbit antibody



Catalog Number: 167591

## Product name

Anti-Ly-6A/E Rabbit antibody

## Specificity

Mouse

## Antibody description

Ly-6A/E Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of mouse Ly-6A/E

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/50

## Validations



Western blot detection of Ly-6A/E in HeLa,A549,MCF-7,HL-60 cell lysates using Ly-6A/E Rabbit mAb(1:1000 diluted).Predicted band size:14kDa.Observed band size:18kDa.