

Anti-TriMethyl-Histone H3 (Lys27) Rabbit antibody



Catalog Number: 167592

Product name

Anti-TriMethyl-Histone H3 (Lys27) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TriMethyl-Histone H3 (Lys27) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic methyl-peptide corresponding to residues surrounding Lys27 of human Histone H3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, ChIP

Validations



Western blot detection of Histone H3 (tri methyl-Lys27) in C6, Jurkat cell lysates using Histone H3 (tri methyl-Lys27) Rabbit mAb(1:1000 diluted). Predicted band size: 15kDa. Observed band size: 15kDa.