

# Anti-P2Y12 Rabbit antibody



Catalog Number: 167595

## Product name

Anti-P2Y12 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

P2Y12 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human P2Y12

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of P2Y12 in A549,HL-60,U251,U87-MG cell lysates using P2Y12 Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:52kDa.