

Anti-Acetyl-Histone H3 (Lys4) Rabbit antibody



Catalog Number: 167606

Product name

Anti-Acetyl-Histone H3 (Lys4) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB, ICC/IF

Antibody description

Acetyl-Histone H3 (Lys4) Rabbit monoclonal antibody

Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

Preparation

Antigen: A synthesized peptide derived from human Histone H3 (acetyl K4)

Validations

Formulation


Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of Acetyl-Histone H3 (Lys4) in A549,HL-60 cell lysates using Acetyl-Histone H3 (Lys4) Rabbit mAb(1:1000 diluted).Observed band size:15kDa.