

# Anti-GRK6 Rabbit antibody



Catalog Number: 167608

## Product name

Anti-GRK6 Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

GRK6 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human GRK6

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000

ICC/IF: 1/50

## Validations



Western blot detection of GRK6 in A549,HL-60,U2OS cell lysates using GRK6 Rabbit mAb(1:1000 diluted).Predicted band size:66kDa.Observed band size:66kDa.