

# Anti-PELP1 Rabbit antibody



Catalog Number: 167610

## Product name

Anti-PELP1 Rabbit antibody

## Specificity

Human

## Antibody description

PELP1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human PELP1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/1000

ICC/IF: 1/20

FC: 1/50

## Validations



Western blot detection of PELP1 in HeLa,A549,HL-60 cell lysates using PELP1 Rabbit mAb(1:1000 diluted).Predicted band size:120kDa.Observed band size:160kDa.