

Anti-DAB2 Rabbit antibody



Catalog Number: 167611

Product name

Anti-DAB2 Rabbit antibody

Specificity

Human

Antibody description

DAB2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human DAB2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of DAB2 in HeLa, A549, HL-60, U2OS, C6 cell lysates using DAB2 Rabbit mAb(1:1000 diluted). Predicted band size: 83kDa. Observed band size: 96kDa.