

Anti-BarX1 Rabbit antibody



Catalog Number: 167619

Product name

Anti-BarX1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BarX1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human BarX1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of BarX1 in HeLa,CHO-K1,C6,Ramos cell lysates using BarX1 Rabbit mAb(1:1000 diluted).Predicted band size:27kDa.Observed band size:27kDa.