

# Anti-WWOX Rabbit antibody



Catalog Number: 167630

## Product name

Anti-WWOX Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

WWOX Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human WWOX

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/100

FC: 1/50

IP: 1/20

## Validations



Western blot detection of WWOX in HeLa,U2OS,U251MG cell lysates using WWOX Rabbit mAb(1:1000 diluted).Predicted band size:47kDa.Observed band size:47kDa.