

Anti-Glutathione Peroxidase 1 Rabbit antibody



Catalog Number: 167645

Product name

IgG

Anti-Glutathione Peroxidase 1 Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC, IP

Human, Rat

Dilutions

Antibody description

WB: 1/1000

Glutathione Peroxidase 1 Rabbit monoclonal antibody

IHC: 1/50

ICC/IF: 1/20

Preparation

FC: 1/20

Antigen: A synthetic peptide of human Glutathione Peroxidase 1

IP: 1/20

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of Glutathione Peroxidase 1 in A549,HL-60 using Glutathione Peroxidase 1 Rabbit mAb(1:1000 diluted)

