

Anti-TFAP2C/AP2 gamma Rabbit antibody



Catalog Number: 167647

Product name

Anti-TFAP2C/AP2 gamma Rabbit antibody

Specificity

Human

Antibody description

TFAP2C/AP2 gamma Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human AP2 gamma/TFAP2C

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Validations



Western blot detection of AP2 gamma/TFAP2C in HeLa,A549 using AP2 gamma/TFAP2C Rabbit mAb(1:1000 diluted)