

Anti-IQGAP2 Rabbit antibody



Catalog Number: 167652

Product name

Anti-IQGAP2 Rabbit antibody

Specificity

Human, Mouse

Antibody description

IQGAP2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human IQGAP2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/20

Validations



Western blot detection of IQGAP2 in A549,HL-60 using IQGAP2 Rabbit mAb(1:1000 diluted)