

Anti-Entactin/NID Rabbit antibody



Catalog Number: 167666

Product name

Rabbit IgG

Anti-Entactin/NID Rabbit antibody

Applications

Specificity

WB, IHC-P

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

Entactin/NID Rabbit monoclonal antibody

IHC: 1/100

Preparation

Antigen: Recombinant protein of human Entactin/NID

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Entactin/NID in HeLa,A549,HL-60,U251,U87-MG cell lysates using Entactin/NID Rabbit mAb(1:1000 diluted).Predicted band size:136kDa.Observed band size:136kDa.