

Anti-MAP3K4 Rabbit antibody



Catalog Number: 167671

Product name

Anti-MAP3K4 Rabbit antibody

Specificity

Human

Antibody description

MAP3K4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MAP3K4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/20

ICC/IF: 1/50

FC: 1/20

Validations



Western blot detection of MAP3K4 in HL-60,U251 cell lysates using MAP3K4 Rabbit mAb(1:1000 diluted).Predicted band size:182kDa.Observed band size:200kDa.