

# Anti-CD11a Rabbit antibody



Catalog Number: 167672

## Product name

Rabbit IgG

Anti-CD11a Rabbit antibody

## Applications

## Specificity

WB, IHC-P

Mouse

## Dilutions

## Antibody description

WB: 1/1000

CD11a Rabbit monoclonal antibody

IHC: 1/500

## Preparation

Antigen: Recombinant protein of mouse CD11a

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of CD11a in A549,HL-60 cell lysates using CD11a Rabbit mAb(1:1000 diluted).Predicted band size:128kDa.Observed band size:160kDa.