

Anti-Pleckstrin Rabbit antibody



Catalog Number: 167675

Product name

Anti-Pleckstrin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Pleckstrin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Pleckstrin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Pleckstrin in HeLa,HL-60 cell lysates using Pleckstrin Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:40kDa.