

Anti-BAP1 Rabbit antibody



Catalog Number: 167695

Product name

Anti-BAP1 Rabbit antibody

Specificity

Human

Antibody description

BAP1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human BAP1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of BAP1 in HeLa,A549,HL-60,U2OS cell lysates using BAP1 Rabbit mAb(1:1000 diluted).Predicted band size:80kDa.Observed band size:95kDa.