

Anti-TFPI Rabbit antibody



Catalog Number: 167696

Product name

Anti-TFPI Rabbit antibody

Specificity

Human

Antibody description

TFPI Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human TFPI

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000

FC: 1/50

IP: 1/20

Validations



Western blot detection of TFPI in A549,HL-60,U2OS cell lysates using TFPI Rabbit mAb(1:1000 diluted).Predicted band size:35kDa.Observed band size:40-55kDa.