

Anti-Flotillin 1 Rabbit antibody



Catalog Number: 167699

Product name

Anti-Flotillin 1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Flotillin 1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Flotillin 1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Flotillin 1 in HeLa, A549, MCF-7, U251 cell lysates using Flotillin 1 Rabbit mAb (1:1000 diluted). Predicted band size: 47 kDa. Observed band size: 49 kDa.