

Anti-Collagen I Rabbit antibody



Catalog Number: 167705

Product name

Anti-Collagen I Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Collagen I Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Collagen I

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

;ICC/IF: 1/50

;FC: 1/50

;IP: 1/20

Validations



Western blot detection of Collagen I in A549,HL-60,U2OS using Collagen I Rabbit mAb(1:1000 diluted)