

Anti-MRPS18B Rabbit antibody



Catalog Number: 167706

Product name

Rabbit IgG

Anti-MRPS18B Rabbit antibody

Applications

Specificity

WB, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

MRPS18B Rabbit monoclonal antibody

ICC/IF: 1/200

Preparation

IP: 1/20

Antigen: Recombinant protein of human MRPS18B

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of MRPS18B in A549,HL-60,U251 cell lysates using MRPS18B Rabbit mAb(1:1000 diluted).Predicted band size:29kDa.Observed band size:29kDa.