

Anti-BCRP/ABCG2 Rabbit antibody



Catalog Number: 167711

Product name

Anti-BCRP/ABCG2 Rabbit antibody

Specificity

Human

Antibody description

BCRP/ABCG2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human BCRP/ABCG2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of BCRP/ABCG2 in HepG2,U251,HL-60 using BCRP/ABCG2 Rabbit mAb(1:1000 diluted)