

Anti-Noxa Rabbit antibody



Catalog Number: 167719

Product name

Anti-Noxa Rabbit antibody

Specificity

Human

Antibody description

Noxa Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Noxa

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Noxa in HeLa,A549,HL-60,U251 cell lysates using Noxa Rabbit mAb(1:1000 diluted).Predicted band size:6kDa.Observed band size:10kDa.