

Anti-NDUFB1 Rabbit antibody



Catalog Number: 167737

Product name

Anti-NDUFB1 Rabbit antibody

Specificity

Human

Antibody description

NDUFB1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFB1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000-1/2000

IHC: 1/100

Validations



Western blot detection of NDUFB1 in A549,HL-60,U2OS cell lysates using NDUFB1 Rabbit mAb(1:1000 diluted).Predicted band size:7kDa.Observed band size:7kDa.