

Anti-PACSIN3 Rabbit antibody



Catalog Number: 167738

Product name

Anti-PACSIN3 Rabbit antibody

Specificity

Human

Antibody description

PACSIN3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human PACSIN3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000-1/5000


IHC: 1/20

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of PACSIN3 in HeLa,A549 cell lysates using PACSIN3 Rabbit mAb(1:1000 diluted).Predicted band size:49kDa.Observed band size:49kDa.