

Anti-HJURP Rabbit antibody



Catalog Number: 167758

Product name

Anti-HJURP Rabbit antibody

Specificity

Human

Antibody description

HJURP Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human HJURP

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/20

Validations



Western blot detection of HJURP in HeLa,CHO-K1, using HJURP Rabbit mAb(1:1000 diluted)