

Anti-Tissue Factor Rabbit antibody



Catalog Number: 167773

Product name

Anti-Tissue Factor Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Tissue Factor Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Tissue Factor

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/50

ICC/IF: 1/50-1/200

Validations



Western blot detection of Tissue Factor in HeLa,A549,MCF-7,U251,HL-60 cell lysates using Tissue Factor Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:45kDa.