

# Anti-Caspase-7 Rabbit antibody



Catalog Number: 167776

## Product name

Rabbit IgG

Anti-Caspase-7 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Caspase-7 Rabbit monoclonal antibody

IHC: 1/200

## Preparation

FC: 1/20

Antigen: Recombinant protein of human Caspase-7

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Caspase-7 in HeLa,A549,MCF-7,U251,HL-60 cell lysates using Caspase-7 Rabbit mAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.