

# Anti-phospho-TrkB (Tyr705) Rabbit antibody



Catalog Number: 167778

## Product name

Monoclonal

Anti-phospho-TrkB (Tyr705) Rabbit antibody

## Ig Type

## Specificity

IgG

Human, Mouse, Rat

## Applications

## Antibody description

WB

phospho-TrkB (Tyr705) Rabbit monoclonal antibody


## Dilutions

WB: 1/1000

## Preparation

## Validations

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr705 of human TrkB



## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Western blot detection of phospho-TrkB (Tyr705) in HeLa,A549,C6 cell lysates using phospho-TrkB (Tyr705) Rabbit mAb(1:1000 diluted).Predicted band size:92kDa.Observed band size:140,90kDa.