

Anti-Ku70 Rabbit antibody



Catalog Number: 167782

Product name

Anti-Ku70 Rabbit antibody

Specificity

Human

Antibody description

Ku70 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Ku70

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/50

Validations



Western blot detection of Ku70 in HeLa, A549, MCF-7, U251 cell lysates using Ku70 Rabbit mAb(1:1000 diluted). Predicted band size: 70kDa. Observed band size: 70kDa.