

Anti-SQSTM1/p62 Rabbit antibody



Catalog Number: 167788

Product name

Anti-SQSTM1/p62 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SQSTM1/p62 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human SQSTM1/p62

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of SQSTM1/p62 in A549,HL-60,U2OS cell lysates using SQSTM1/p62 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:62kDa.