

# Anti-UBE3A Rabbit antibody



Catalog Number: 167789

## Product name

Anti-UBE3A Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

UBE3A Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human UBE3A

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

## Validations



Western blot detection of UBE3A in A549,HL-60,U2OS cell lysates using UBE3A Rabbit mAb(1:1000 diluted).Predicted band size:101kDa.Observed band size:101kDa.