

# Anti-SV40 T-antigen Rabbit antibody



Catalog Number: 167790

## Product name

Anti-SV40 T-antigen Rabbit antibody

## Specificity

ALL

## Antibody description

SV40 T-antigen Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of SV40 T-antigen

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/50

IP: 1/20

## Validations



Western blot detection of SV40 T-antigen in A549,HL-60,U251 cell lysates using SV40 T-antigen Rabbit mAb(1:1000 diluted).Predicted band size:82kDa.Observed band size:90kDa.