

Anti-NUP50 Rabbit antibody



Catalog Number: 167792

Product name

Anti-NUP50 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NUP50 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NUP50

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/50

FC: 1/20-1/100

Validations



Western blot detection of NUP50 in HeLa,C6 cell lysates using NUP50 Rabbit mAb(1:1000 diluted).Predicted band size:50kDa.Observed band size:50kDa.