

Anti-NUP62CL Rabbit antibody



Catalog Number: 167793

Product name

Anti-NUP62CL Rabbit antibody

Specificity

Human

Antibody description

NUP62CL Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NUP62CL

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of NUP62CL in HeLa, Ramos cell lysates using NUP62CL Rabbit mAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 21kDa.