

# Anti-MTAP Rabbit antibody



Catalog Number: 167798

## Product name

Anti-MTAP Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MTAP Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human MTAP

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

## Validations



Western blot detection of MTAP in A549, MCF-7, U251, HL-60 cell lysates using MTAP Rabbit mAb(1:1000 diluted). Predicted band size: 31kDa. Observed band size: 31kDa.