

Anti-Bora Rabbit antibody



Catalog Number: 167799

Product name

Anti-Bora Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bora Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Bora

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of Bora in HL-60,C6 cell lysates using Bora Rabbit mAb(1:1000 diluted).Predicted band size:61kDa.Observed band size:61kDa.