

Anti-Sortilin/NT3 Rabbit antibody



Catalog Number: 167814

Product name

Anti-Sortilin/NT3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Sortilin/NT3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Sortilin/NT3

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/500

ICC/IF: 1/20

FC: 1/50

Validations



Western blot detection of Sortilin/NT3 in HeLa, A549, MCF-7 cell lysates using Sortilin/NT3 Rabbit mAb(1:1000 diluted). Predicted band size: 92kDa. Observed band size: 92kDa.